

五、课题组管理 PI GROUP MANAGEMENT

测试课题组
P.I.: 仪器管理员

Edit

1. Find PI's Laboratory page, check the "Inactive user" label;

Laboratory Members 预约审核 Access History Equipment Log Equipment Sign-Up Billing Testing Service Accounts

Current User 1 Inactive User

Search	Name	Contact Info	From - To	Address	
	高胜寒 Inactive	123123123 5df787c877711	unknown - 12/16/2019		2. Edit Follow

2. Click "Edit" for activating the user;

五、课题组管理 PI GROUP MANAGEMENT

Get the page of user's profile:

People 高胜寒 Edit

General Picture Roles Account Equipment Sign-Up Laboratory Notification 负责门禁

不审核
Send a message

1. Personal information

Name 高胜寒 *required
Gender Male
Title PhD Candidate *required
Group All +
ID Number 5df787c877711
Major
Affiliation
From - To 2020/02/05 - 2019/12/16

2. Contact Info

Email 5df787c877711 *required
Phone 123123123 *required
Address

IC Card No.

Activated Yes No
Unremovable Yes No

Update Delete

3. Choose “Yes” for “Activated” option;

4. Click “Update” for saving the setting.